南模生物TM
S H A N G H A I
M O D E L O R G A N I S M S

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
Wang MS, Huo YX, Li Y, et al. Comparative population genomics reveals genetic basis underlying body size of domestic chickens. **J Mol Cell Biol.** 2016 Dec;8(6):542-552.

Feng N, Chen H, Fu S, et al. HIF-1 α and HIF-2 α induced angiogenesis in gastrointestinal vascular malformation and reversed by thalidomide. **Sci Rep.** 2016 Jun 1;6:27280.

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Wang MS, Zhang RW, Su LY, et al. Positive selection rather than relaxation of functional constraint drives the evolution of vision during chicken domestication. **Cell Res.** 2016 May;26(5):556-73.

Gao X, Yuan YY, Lin QF, et al. Mutation of IFNLR1, an interferon lambda receptor 1, is associated with autosomal-dominant non-syndromic hearing loss. **J Med Genet.** 2018 Feb 16.

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Morpholino Knockdown 靶点位置

图 3. Morpholino Knockdown 靶点位置示意图。 The zebrafish gene-X was targeted by specific morpholino antisense strategies to prevent either the translation of the zebrafish gene (ATG-MO) or proper splicing of exon 3 (E3I3-MO).

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
图 4. Panels **A** through **H** show lateral views of control MO injected zebrafish embryos (Panel **A** and Panel **E**) and embryos injected with gene-X morpholino oligonucleotides (MO) (Panel **B** through **C**, Panel **F** through **G**), gene-X-e3i3-MO plus nonmutant zebrafish gene-X (Panel **D** and Panel **H**). Coinjection of nonmutant zebrafish gene-X mRNA rescued U-shaped somites (red arrow) and curved body axis (blue dotted line) in gene-X morphants at 52 hpf. The bar graph in Panel **I** shows the percentage of embryos with development defects. . Panels **J** Effectiveness of gene-X knockdown was confirmed by RT-PCR and sanger sequencing. hpf, hours post fertilization.

Morpholino Knockdown 信号通路机制分析


图 5. 基因 Knockdown 对 Hedgehog 信号通路影响。Endogenous shha, ptch1, ptch2, sufu, gli1, gli2a, gli2b, and gli3 in wild-type control and geneX morphants assessed by qRT-PCR (n = 100 individual embryos). GeneX fine-tunes Hedgehog patterning activity maybe through a novel regulatory feedback loop.

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
图 6. Morpholino knock down of geneX induces potent CNS-specific apoptosis. Wild-type control embryos and embryos injected with geneX-MO were stained with acridine orange (AO) at 32hpf. Apoptotic cells are visible as black spots, and less bright homogenous black staining is unspecific background staining. (A–B) Uninjected wild-type control zebrafish exhibited few or no apoptotic cells in CNS (central nervous system). In contrast, significantly increased staining was observed throughout the CNS in zebrafish injected with geneX-MO (C–D). The blue boxed regions are shown at higher magnification in the right panels. A–D: lateral view, anterior, left.

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
图 7. Morpholino knock down of geneX induces apoptosis. geneX knock down induces potent ototoxicity in zebrafish. Wild-type control embryos and embryos injected with geneX-MO were stained with the mitochondrial potentiometric dye DASPEI at 6-dpf. Hair cells stereotypically located on the lateral line were stained as green dots (white arrow). Uninjected wild-type control zebrafish exhibited normal hair cell number. In contrast, significantly decreased hair cell staining was observed in zebrafish injected with geneX-MO. Fluorescent DASPEI images were inverted for particle analysis. The fluorescence particle signal was quantified using morphometric analysis. dpf, days post fertilization.

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
图 8. GeneX knock down inhibits the trunk angiogenesis in zebrafish. (A–D) Representative fluorescent images of zebrafish embryos at 32h post-fertilization (hpf). (C–E) Compared with wild-type control, embryos injected with geneX-MO present a lower number of incomplete ISVs and only occasional sprouts (asterisk) of dorsal aorta. The boxed regions are shown at higher magnification in the right panels. DLAV, dorsal longitudinal anastomotic vessels; ISV, intersegmental vessels; DA, dorsal aorta; PCV, posterior cardinal vein.

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Tu X, Deng Y, Chen J, et al. Screening study on the anti-angiogenic effects of Traditional Chinese Medicine - Part I: heat-clearing and detoxifying TCM. **J Ethnopharmacol.** 2016 Sep 8.

Guan YY, Liu HJ, Luan X, et al. Raddeanin A, a triterpenoid saponin isolated from Anemone raddeana, suppresses the angiogenesis and growth of human colorectal tumor by inhibiting VEGFR2 signaling. **Phytomedicine.** 2015 Jan 15;22(1):103-10.

Zhang J, Gao B, Zhang W, et al. Monitoring antiangiogenesis of bevacizumab in zebrafish. **Drug Des Devel Ther.** 2018 Aug 6;12:2423-2430.

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Sun Y, Sheng Z, Ma C, et al. Combining genomic and network characteristics for extended capability in predicting synergistic drugs for cancer. **Nat Commun.** 2015 Sep 28;6:8481.

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
Zhang S, Zhang Q, Zhang D, et al. Anti-osteoporosis activity of a novel *Achyranthes bidentata* polysaccharide via stimulating bone formation. **Carbohydr Polym.** 2018 Mar 15;184:288-298.

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
图 18. GeneX knockdown induced a reduction of the distance that embryos covered in 5 min compared to control.

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MO knock-down	Wang MS, Zhang RW, Su LY, et al. Positive selection rather than relaxation of functional constraint drives the evolution of vision during chicken domestication. Cell Res. 2016 May;26(5):556-73. Gao X, Yuan YY, Lin QF, et al. Mutation of IFNLR1, an interferon lambda receptor 1, is associated with autosomal-dominant non-syndromic hearing loss. J Med Genet. 2018 Feb 16.
Over-expression	Du Z, Ma HL, Zhang ZY, et al. Transgenic Expression of A Venous Malformation Related Mutation, TIE2-R849W, Significantly Induces Multiple Malformations of Zebrafish. Int J Med Sci. 2018 Feb 12;15(4):385-394. Wang MS, Huo YX, Li Y, et al. Comparative population genomics reveals genetic basis underlying body size of domestic chickens. J Mol Cell Biol. 2016 Dec;8(6):542-552. Feng N, Chen H, Fu S, et al. HIF-1 α and HIF-2 α induced angiogenesis in gastrointestinal vascular malformation and reversed by thalidomide. Sci Rep. 2016 Jun 1;6:27280.
抗癌药效学评价	Sun Y, Sheng Z, Ma C, et al. Combining genomic and network characteristics for extended capability in predicting synergistic drugs for cancer. Nat Commun. 2015 Sep 28;6:8481.
抗骨质疏松药药效学评价	Zhang S, Zhang Q, Zhang D, et al. Anti-osteoporosis activity of a novel Achyranthes bidentata polysaccharide via stimulating bone formation. Carbohydr Polym. 2018 Mar 15;184:288-298.
抗血管生成药物药效学评价及筛选	Tu X, Deng Y, Chen J, et al. Screening study on the anti-angiogenic effects of Traditional Chinese Medicine - Part I: heat-clearing and detoxicating TCM. J Ethnopharmacol. 2016 Sep 8. Guan YY, Liu HJ, Luan X, et al. Raddeanin A, a triterpenoid saponin isolated from Anemone raddeana, suppresses the angiogenesis and growth of human colorectal tumor by inhibiting VEGFR2 signaling. Phytomedicine. 2015 Jan 15;22(1):103-10. Zhang J, Gao B, Zhang W, et al. Monitoring antiangiogenesis of bevacizumab in zebrafish. Drug Des Devel Ther. 2018 Aug 6;12:2423-2430.
毒理学评价	Wang Q, Liu S, Hu D, et al. Identification of apoptosis and macrophage migration events in paraquat-induced oxidative stress using a zebrafish model. Life Sci. 2016 Jul 15;157:116-24.

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Shanghai ZJ Biomedical Service Alliance

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